

SmartBioSurface® Slides

Where efficient cell adhesion meets easy isolation

Single-Cell Laser Capture Microdissection of Non-Adherent Cells

Single-cell Laser Capture Microdissection (LCM) enables high-resolution molecular profiling but requires **efficient immobilization and release**.

- **Non-adherent cells** (PBMCs, suspension cells, and rare cell populations) are **difficult to immobilize** on standard LCM slides, and cytopsin leads to low recovery and cell damage
- **SmartBioSurface® slides**, thanks to nanostructured TiO₂ coating, **promote cell adhesion**, preserve **morphology**, enable **laser microdissection**, and **downstream molecular analysis**



SmartBioSurface® slides. Left Image: SmartBioSurface® slides in a different format. Right Image: Atomic Force Microscopy (AFM) of the nanocoating (50 nm thickness)

LCM Workflow on SmartBioSurface® Slides






Schematic description of the analytical workflow starting from cell seeding on SmartBioSurface® slides to single-cell genomic analysis.

- **Cell seeding** – Non-adherent living cells are directly pipetted onto the slide and spontaneously adhere to the nanostructured surface
- **Staining** – Samples can be processed with brightfield or immunofluorescence staining for clear cell identification
- **LCM isolation** – Target cells are precisely identified and isolated using laser microdissection
- **Molecular analysis** – Recovered single cells are compatible with WGA, NGS library preparation, and CNV analysis

Easy and Efficient Laser Microdissection

SmartBioSurface® slides enable **clear target identification** and clean **laser-based release** with LCM systems, allowing **efficient single-cell recovery**.

Step 1 Target Identification	Step 2 Post-Dissection	Step 3 Capture Confirmation
		
<i>Target cell identified on SmartBioSurface® slide</i>	<i>Excision Site</i>	<i>Cell recovered in the collection cap</i>

Target cell identification, laser microdissection, and successful cell capture in the collection cap: the white arrow indicates the target undergoing microdissection (Zeiss PALM MicroBeam LCM).

SmartBioSurface® Slides Key Advantages for LCM Applications

- **Ready-to-use** slides with no specialized preparation
- Efficient **adhesion** of non-adherent and rare cells
- Compatible with **brightfield** and **immunofluorescence** staining
- 100% efficiency of single-cell **microdissection**
- Compatible with **downstream molecular analysis**

References

1. Chen S, El-Heliebi A, Schmid J, Kashofer K, Czyż ZT, Polzer BM, Pantel K, Kroneis T, Sedlmayr P. Target Cell Pre-enrichment and Whole Genome Amplification for Single Cell Downstream Characterization. *J Vis Exp.* 2018 May 15;(135):56394. doi: 10.3791/56394. PMID: 29863657; PMCID: PMC6101176
2. Carbone R, Marangi I, Zanardi A, et al. Biocompatibility of cluster-assembled nanostructured TiO₂ with primary and cancer cells. *Biomaterials.* 2006;27(17):3221-3229. doi: 10.1016/j.biomaterials.2006.01.056
3. Zanardi A, Bandiera D, Bertolini F, et al. Miniaturized FISH for screening of onco-hematological malignancies. *Biotechniques.* 2010;49(1):497-504. doi:10.2144/000113445
4. Clinical trial NCT05942066
5. Visci G, Tolomeo D, Lonoce A, et al. A novel method for the isolation of single cells mimicking circulating tumour cells adhered on Smart Bio Surface slides by Laser Capture Microdissection. *PLoS One.* 2024;19(3): e0297739. Published 2024 Mar 8. doi: 10.1371/journal.pone.0297739

Tethis S.p.A.

Via Francesco Olgiati 5, 20143 Milan, Italy

E-mail: contact@tethis-lab.com

Visit us

www.tethis-lab.com